



Nexus screens included with the Audio-Visual Library

www.growing.com/mindware

Neurofeedback Dual Monitor Master



Client-side video feedback (NF or BF)



Client-side DVD feedback (NF or BF)

Neurofeedback Single Monitor Video Left-sided



Neurofeedback Single Monitor Video Right-sided



Neurofeedback Single Monitor DVD Left-sided



Neurofeedback Single Monitor DVD Right-sided



Biofeedback Single Monitor Video Left-sided



Biofeedback Single Monitor Video Right-sided



Biofeedback Single Monitor DVD Left-sided



Biofeedback Single Monitor DVD Right-sided



Biofeedback Dual Monitor Master



Client-side video feedback (NF or BF)



Client-side DVD feedback (NF or BF)

Dual Monitor Screens for Biofeedback and Neurofeedback

There are two "Dual-Monitor Master" screens, one for biofeedback and one for neurofeedback. Each can be used with either "Client-side video feedback" or "Client-side DVD feedback".

Single Monitor Screens for Biofeedback and Neurofeedback

Note that the screens on the left are similar to those on the right. Choose between them based on where your client is seated. The single monitor neurofeedback screens come in two versions, one for video file feedback (row one above) and one for DVD feedback (row two). The same is true for the single monitor biofeedback module shown in rows 3 and 4.