


Low-Frequency Neurofeedback Module for Nexus

For 1 Bipolar Channel - Includes single screen version and dual-screen version



Fully programmable reward module. "Hot-Key" adjustable frequency.

Raw Electrocardiograph

HRV Power Spectrum

Sensor "DC Offset"

Temperature

Galvanic skin response

Heart rate

Heart rate variability

HRV - low frequency %

Raw EEG

Spectral Display

10 programmable inhibit bins

% time above threshold

Includes video sampler to use for feedback on single or dual screen systems.

You can use any commercial DVD or your own custom AVI file for visual feedback.

Comes with audio sampler for feedback. You can also use your own mp3 and wav files.

Comes with : Instructions; Screens, channel set, video and audio sampler files for feedback, user's manual and protocol guide.

Registered users are entitled to free technical support. Upgrades are always free.

Custom Nexus applications can be handcrafted to your specifications. For more information call Douglas Dailey at (408) 984-3333 or email: mind@growing.com